
# **bone**scalpel<sup>®</sup> 超声波骨

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## BoneScalpel Technology Comparison

	Power			Manual	
	BoneScalpel	High Speed Drill	Micro Saw	Osteotome	Kerrison Punch
Cutting frequency	22,500 strokes/sec (1,300,000 strokes/ min)	~ 1,300 rotations/sec (80,000 rpm)	~ 300 cps (20,000 cpm)	single stroke	single punch
Handpiece configuration	straight / angled	straight / angled	straight / angled	straight	angled
Cutting mode	longitudinal	rotational	Longitudinal transvers	longitudinal	jawlike rotation
Tip types	multiple osteotome + shaver + drill (coming)	multiple burr + drill	Single sagittal oscillating saws	single osteotomes	single punches
Tip cooling	direct to active edge	indirect / ancilliary	indirect / ancilliary	n/a	n/a
Precision	very high	high	high	high	high
Minimal kerf size	0.5 mm	2 mm	0.3 mm	2-3 mm	2-3 mm
Axial tip displacement	very low	medium	medium	very high	n/a
Lateral tip walking / straying	very low	very high (walking/skyving)	very high (bilateral straying)	n/a	n/a
Tip inertia	very low	very high	high	very high	n/a
Tip stop/start/stop	near instantenous	delayed	delayed	cannot be stopped after strike	n/a
Impact force	very low	low	medium	very high	high
Risk of tip slippage	very low	very high risk of tissue injury	low	very high	none
En bloc bone dissection for harvest	very good	poor	moderate	very good	n/a
Frontal bone dissection	very good	very good	poor	very good	n/a
Lateral bone dissection	Limited with blade, very good with shaver	very good	good	n/a	moderate
Fine bone ablation	very good	very good	n/a	n/a	n/a
Bone debris	minimal not visible & flushed out	very high	minimal	n/a	very high bite chunks - require assistance with every actuation
Tissue selectivity	very high	n/a	n/a	n/a	n/a
Soft tissue interaction	gentle / atraumatic (reduced on shaving surfaces)	very aggressive	very aggressive	aggressive	aggressive
Tissue grabbing/ tearing	none	very aggressive risk of tissue injury	aggressive on teeth	low	low to none
Ability to use cottonoids	yes	no	no	n/a	n/a
Reduced bone bleeding	generally observable significant in advanced osetotomies	n/a	n/a	n/a	n/a
Thermal consideration	control irrigation & exposure time	control ancilliary irrigation	control ancilliary irrigation	no	no
Repetitive hand fatigue	low	low	low	low	very high

**Table 1.** Comparison of BoneScalpel to other bone cutting tools.

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## References

1. Pakzaban P. BoneScalpel™ Ultrasonic Bone Cutting System: Neurosurgical and Orthopaedic Applications, Aesculap Power Systems, [www.aesculapusa.com](http://www.aesculapusa.com), 2012
2. Sanborn MR, Balzer J, Gerszten PC, Karausky P, Cheng BC, Welch WC. Safety and efficacy of a novel ultrasonic osteotome device in an ovine model. *J Clin Neurosci*. 2011 Nov;18(11):1528-33. Epub 2011 Sep
3. Recinos VR, Ahn E, Carson B, Jallo G. Technical Note: A Novel Bone-Cutting Instrument, the BoneScalpel™, May be Useful in Performing Osteoplastic Laminoplasty. AANS/CNS Section on Pediatric Neurological Surgery Annual Meeting, Boston, MA, December 2009.
4. Parker SL, Kretzer RM, Recinos PF, Jallo GI, Recinos VR. Use of a Novel Ultrasonic Bone Scalpel for Osteoplastic Laminoplasty in the Resection of Intradural Spinal Cord Pathology. 27th Annual Meeting of the AANS/CNS Section on Disorders and Peripheral Nerves, Phoenix, AZ, March 2011
5. Parker SL, Kretzer RM, Recinos PF, Molina CA, Wolinsky JP, Jallo GI, Recinos VR. Ultrasonic Bone Scalpel for Osteoplastic Laminoplasty in the Resection of Intradural Spinal Pathology: Case Series and Technical Note. *Neurosurgery*. 2012 Dec 20.
6. Nickele C, Hanna A, Baskaya MK. Osteotomy for Laminoplasty without Soft Tissue Penetration, Performed Using a Harmonic Bone Scalpel: Instrumentation and Technique. *J Neurol Surg A Cent Eur Neurosurg*. 2013 Mar 15.
7. Lieberman IH, Hu X. Use of Ultrasonic BoneScalpel in Spine Surgeries: Experience from the First 58 Patients. 19th International Meeting on Advanced Spine Techniques (IMAST), Istanbul, Turkey, July, 2012
8. Hu X, Ohnmeiss DD, Lieberman IH. Use of an ultrasonic osteotome device in spine surgery: experience from the first 128 patients. *Eur Spine J*. 2013 Apr 16.
9. Suchomel P, Hradil J. [Minimally invasive cervical elastic laminoplasty - principles and surgical technique]. *Acta Chir Orthop Traumatol Cech*. 2011;78(5):437-41. Czech.
10. Chalchana KL, Jallo GI, Dorafshar AH, Ahn ES. Novel use of an ultrasonic bone-cutting device for endoscopic-assisted craniostylosis surgery. *Childs Nerv Syst*. 2013 Feb 6: DOI 10.1007/s00381-013-2043-6. [Epub ahead of print].
11. Gilles R, Couvreur T, Dammous S. Ultrasonic orthognathic surgery: enhancements to established osteotomies. *Int J Oral Maxillofac Surg*. 2013 Jan 8.
12. Wick CC, Rezaee RP, Zender CA. Piezoelectric BoneScalpel osteotomies in osteocutaneous free flaps. *Laryngoscope*. 2013 Mar;123(3):618-21. Epub 2012 Sep 24.
13. Dammous S, Gilles R. Use of Ultrasonic Straight Blade in Orthognathic Surgery: Review of 75 Patients. *Int J Oral Maxillofac Surg*. 2011 40 (10): 1082.
14. Pakzaban P. BoneScalpel™ in Spine Surgery. Surgical Technique and Operative Approaches, Aesculap Power Systems, [www.aesculapusa.com](http://www.aesculapusa.com), 2011

1938 NEW HWY, FARMINGDALE, N.Y. 11735 | +1.631.694.9555 +1.631.694.3285 FAX  
MISONIX, INC. | NASDAQ SYMBOL. MSON | MISONIX.COM