Which **Products** for Which **Targets**

Products	Approaches	Procedures	Notes & Comments
10mm & 20mm Std. (MXB-10 & MXB-20)	Open Mini-Open	Laminectomy Laminotomy Foraminotomy Microdiscectomy TLIF / PLIF Facetectomy	Standard Probe Workhorse – Most Efficient 10mm = Faster Cutting Selective Mini-Open not all
10mm & 20mm LC (MXB-10LC & MXB-20LC)	Tubular Retractor Rigid Tube	Laminectomy Laminotomy Foraminotomy Microdiscectomy TLIF / PLIF Facetectomy	Longer Working Length Increased Visibility Consideration for exposed tissue in tubular retractors Experienced Users – Less Control Recommended Tube(s): 18-22mm Diameter +8cm in Length Rigid Tubes
10mm MIS 20mm MIS Micro Hook MIS (MXB-MIS-10, MXB-MIS-20LCRS, MXB-MIS-S1)	Open Mini-Open Tubular Retractor Rigid Tube	Laminectomy Laminotomy Foraminotomy Microdiscectomy TLIF / PLIF Facetectomy	Protection of non-targeted tissue Most Versatile Product Option Allows for Distal Gripping Positions Recommended Tube(s): 22-26mm Diameter 6-8cm in Length Expandable/NuVasive Tubes

